


Comparison of Shodex Reversed Phase Chromatography (RPC) Column Features

ODS columns are the most popular reversed phase columns that are packed with silica-based octadecyl group. Shodex provides not only ODS columns but also polymer-based reversed phase columns with different functional groups. Please use following descriptions about the column features as guidelines to select suitable columns for your application purposes.


Features

ODP2 HP	<ul style="list-style-type: none">• Provides a large theoretical plate number nearly twice as much as generally available polymer-based reversed phase columns do• Offers enhanced retention of high polar substances compared to ODS columns• Suitable for the analysis of small molecules such as pharmaceuticals in the presence of protein matrix• Ideal for LC/MS analysis of high polar compounds• Fulfills USP L39 requirements
ODP-50	<ul style="list-style-type: none">• Relatively large pore size is suitable for the analysis of amino acids, peptides, and proteins
C8P-50	<ul style="list-style-type: none">• Usable in a wide pH range from pH 2 to 13
C4P-50	<ul style="list-style-type: none">• Usable in 100 % water and buffer solution• Best used for the analysis of basic substances• ODP-50 fulfills USP L67 requirements
ODP-40	<ul style="list-style-type: none">• Higher performance type of ODP-50 series• Fulfills USP L67 requirements
RP18-415	<ul style="list-style-type: none">• Large pore size is suitable for the analysis of proteins and peptides• Fulfills USP L21 requirements
DS-613	<ul style="list-style-type: none">• Suitable for reversed phase analysis of highly hydrophilic substances that are not well retained by ODS columns
DS-413	<ul style="list-style-type: none">• Fulfill USP L21 requirements
DE-613	<ul style="list-style-type: none">• General purpose polymer-based column having similar polarity as ODS columns
DE-413	<ul style="list-style-type: none">• Wide working pH range (from pH 2 to 12), usable in 100 % water and buffer solutions
DE-213	<ul style="list-style-type: none">• Fulfill USP L71 requirements
DM-614	<ul style="list-style-type: none">• Suitable for the analysis of amino acids and water-soluble vitamins• Fulfills USP L39 requirements
NN-814	<ul style="list-style-type: none">• The packing material modified with sulfo groups supports multimode (reversed phase and cation exchange) analysis
NN-614	<ul style="list-style-type: none">• Ideal for the analysis of complex samples containing neutral and ionic substances
NN-414	<ul style="list-style-type: none">• Ideal for the analysis of complex samples containing neutral and ionic substances
JJ-50	<ul style="list-style-type: none">• The packing material is modified with trace amounts of quaternary ammonium groups, and supports multimode (reversed phase and anion exchange) analysis• Ideal for analysis of complex samples containing neutral and ionic substances
C18	<ul style="list-style-type: none">• Fully end capped ODS column available at very reasonable price• Fulfills USP L1 requirements
C18M	<ul style="list-style-type: none">• Monomeric type ODS column, fully end capped high purity silica (99.99 % or higher)• Fulfills USP L1 requirements
C18P	<ul style="list-style-type: none">• Polymeric type ODS column, fully end capped high purity silica (99.99 % or higher)• Excellent acid tolerance• Advantageous for separating planar and nonplanar compounds from each other• Fulfills USP L1 requirements
New C18U	<ul style="list-style-type: none">• UHPLC column (Maximum pressure: 100 MPa)• Achieves high performance analysis with sub-2 μm particles• Organic/inorganic silica hybrid particles provide excellent resolution and mechanical stability and improved alkali durability (from pH 1 to 12)• Usable in 100 % water and buffer solution• Fulfills USP L1 requirements
5C8	<ul style="list-style-type: none">• Use when the retention capacity of C18 is too strong• Rapid mass transfer and fast equilibration allow its use as an ion-pair chromatography• Fulfills USP L7 requirements
5CN	<ul style="list-style-type: none">• Utilizes reversed phase interaction and π-electron interaction to separate regioisomers, which typically cannot be separated with ODS or C8 columns• Fulfills USP L10 requirements
5NPE	<ul style="list-style-type: none">• Utilizes several types of interactions based on π-electrons to separate structural isomers

The interrelation between hydrophobicity and retention, and the interrelation between steric selectivity and retention were compared among Shodex columns for reversed phase chromatography. The retention factor (k') of amylbenzene was used as the retention, the separation factor (α) between n-butylbenzene and n-amino benzene was used as the hydrophobicity. The separation factor between o-terphenyl and triphenylene was used as the steric recognition. Larger separation factor means higher hydrophobicity and higher steric selectivity.




Hydrophobicity differences among Shodex RPCs




Column size : 4.6 mm I.D. x 150 mm each
 Eluent : H₂O/CH₃OH=20/80
 Flow rate : 1.0 mL/min
 Detector : UV (254 nm)
 Column temp. : 40 °C

Steric selectivity differences among Shodex RPCs




Column size : 4.6 mm I.D. x 150 mm each
 Eluent : H₂O/CH₃OH=20/80
 Flow rate : 1.0 mL/min
 Detector : UV (254 nm)
 Column temp. : 40 °C

Comparison of different functional groups on the separation of alkylalcohols



Effects of steric selectivity differences



Polymer-based Reversed Phase Chromatography Columns (RSpak)

Please refer to “Comparison of Shodex Reversed Phase Chromatography (RPC) Column Features” on page 6 and 7 for features.

● Standard columns

Product Code	Product Name	Plate Number (TP/column)	Base Material	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7009000	RSpak RP18-415	≥ 5,000	Styrene divinylbenzene copolymer	6	450	4.6 x 150	H ₂ O/CH ₃ CN=5/95
F6709558	RSpak RP18-G	(guard column)	Styrene divinylbenzene copolymer	6	–	4.6 x 10	H ₂ O/CH ₃ CN/THF=40/30/30
F7001001	RSpak DS-613	≥ 6,500	Styrene divinylbenzene copolymer	6	200	6.0 x 150	H ₂ O/CH ₃ CN/THF=30/40/30
F7001012	RSpak DS-413	≥ 11,000	Styrene divinylbenzene copolymer	3.5	200	4.6 x 150	H ₂ O/CH ₃ CN/THF=40/30/30
F6700140	RSpak DS-G	(guard column)	Styrene divinylbenzene copolymer	10	–	4.6 x 10	H ₂ O/CH ₃ CN/THF=30/40/30
F7001004	RSpak DE-613	≥ 7,000	Polymethacrylate	6	25	6.0 x 150	H ₂ O
F7001005	RSpak DE-413	≥ 11,000	Polymethacrylate	4	25	4.6 x 150	H ₂ O/CH ₃ CN=50/50
F7009030	RSpak DE-413L	≥ 17,000	Polymethacrylate	4	25	4.6 x 250	H ₂ O/CH ₃ CN=50/50
F6700150	RSpak DE-G 4A	(guard column)	Polymethacrylate	10	–	4.6 x 10	H ₂ O
F7001002	RSpak DM-614	≥ 4,500	Polyhydroxymethacrylate	10	200	6.0 x 150	5 mM H ₃ PO ₄ aq.
F6700160	RSpak DM-G 4A	(guard column)	Polyhydroxymethacrylate	12	–	4.6 x 10	5 mM H ₃ PO ₄ aq.

● Semi-micro columns

Product Code	Product Name	Plate Number (TP/column)	Base Material	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7001007	RSpak DE-213	≥ 8,000	Polymethacrylate	4	25	2.0 x 150	H ₂ O/CH ₃ CN=50/50
F6700151	RSpak DE-G 2A	(guard column)	Polymethacrylate	6	–	2.0 x 10	H ₂ O/CH ₃ CN=50/50


● Standard columns

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Base Material	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7008140	RSpak NN-814	≥ 9,000	Sulfo	Polyhydroxymethacrylate	10	200	8.0 x 250	0.1 M Sodium phosphate buffer (pH3.0)
F7008150	RSpak NN-614	≥ 4,000	Sulfo	Polyhydroxymethacrylate	10	200	6.0 x 150	0.1 M Sodium phosphate buffer (pH 3.0)
F6700510	RSpak NN-G	(guard column)	Sulfo	Polyhydroxymethacrylate	10	-	6.0 x 50	0.1 M Sodium phosphate buffer (pH 3.0)
F7008160	RSpak NN-414	≥ 6,000	Sulfo	Polyhydroxymethacrylate	10	200	4.6 x 150	0.1 M Sodium phosphate buffer (pH 3.0)
F7008240	RSpak JJ-50 4D	≥ 4,500	Quaternary ammonium	Polyvinyl alcohol	5	100	4.6 x 150	H ₂ O/CH ₃ CN=40/60

● Semi-micro columns

Product Code	Product Name	Plate Number (TP/column)	Functional Group	Base Material	Particle Size (µm)	Pore Size (Å)	Column Size (mm) I.D. x Length	Shipping Solvent
F7008220	RSpak JJ-50 2D	≥ 3,500	Quaternary ammonium	Polyvinyl alcohol	5	100	2.0 x 150	H ₂ O/CH ₃ CN=40/60


Separation and recovery rate of standard proteins



Recovery (%)	
1. Ribonuclease A	93
2. Insulin	98
3. Cytochrome C	100
4. Lysozyme	100
5. BSA	98
6. Myoglobin	108
7. Ovalbumin	-

Column : Shodex RSpak RP18-415
Eluent : (A); 0.1 % TFA aq./CH₃CN=99/1
 (B); 0.1 % TFA aq./CH₃CN=5/95
 Linear gradient; (B %) 20 % to 60 %, 20 min
Flow rate : 1.0 mL/min
Detector : UV (220 nm)
Column temp. : Room temp.


Alkylbenzenes



Sample : 5 µL
 1. m-Cresol 0.1 %
 2. 2,4-Xylenol 0.1 %
 3. Benzene 0.5 %
 4. Toluene 0.5 %
 5. Ethylbenzene 0.5 %
 6. n-Propylbenzene 0.5 %

Column : Shodex RSpak DS-613
Eluent : H₂O/CH₃CN/THF=30/40/30
Flow rate : 1.0 mL/min
Detector : UV (254 nm)
Column temp. : 40 °C


Fatty acid methyl esters



Sample : 0.2 % each, 20 µL
 1. Methyl linoleate
 2. Methyl palmitate
 3. Methyl oleate
 4. Methyl stearate

Column : Shodex RSpak DS-413
Eluent : H₂O/CH₃CN/THF=25/45/30
Flow rate : 1.0 mL/min
Detector : RI
Column temp. : 40 °C


Organic acids



Sample : 5 µL
 1. Glyoxylic acid 1.78 mg/mL
 2. Tartaric acid 1.95 mg/mL
 3. Malic acid 2.06 mg/mL
 4. Lactic acid 2 µL/mL
 5. Malonic acid 1.95 mg/mL
 6. Acetic acid 2 µL/mL
 7. Succinic acid 2.05 mg/mL
 8. Levulinic acid 1.95 mg/mL
 9. Propionic acid 2 µL/mL

Column : Shodex RSpak DE-413
Eluent : 10 mM H₃PO₄ aq.
Flow rate : 1.0 mL/min
Detector : RI
Column temp. : 50 °C


Food additives (Preservatives)



Sample : 10 µL
 1. Saccharin sodium 0.005 %
 2. p-Hydroxybenzoic acid 0.005 %
 3. Sorbic acid 0.02 %
 4. Benzoic acid 0.02 %
 5. Methyl paraben 0.01 %
 6. Dehydroacetic acid 0.01 %
 7. Ethyl paraben 0.02 %
 8. Propyl paraben 0.02 %

Column : Shodex RSpak DE-413
Eluent : 50 mM KH₂PO₄ + 0.1 % H₃PO₄ aq./CH₃CN
 =65/35
Flow rate : 1.0 mL/min
Detector : UV (210 nm)
Column temp. : 40 °C

Vitamins



Sample : 5 mg/L each (in 2.5 mM H₃PO₄ aq.), 4 µL
 1. Vitamin C
 2. Orotic acid
 3. Vitamin B₆
 4. Vitamin B₁
 5. Nicotinamide
 6. Caffeine

Column : Shodex RSpak DM-614
Eluent : 0.055 M Na₂HPO₄ + 0.045 M KH₂PO₄ aq.
Flow rate : 1.0 mL/min
Detector : UV (254 nm)
Column temp. : 30 °C

