

Shodex Super Summer Sale

We're making room in our warehouse! Selected Shodex columns are now available at heavily reduced prices for a limited time.

Sale period: 1 July – 31 August 2026. While stocks last.

These products have older production dates and are beyond the typical recommended shelf-life period. Product performance should therefore be verified by the customer before use, especially for critical applications.

However, the column chemistry of many Shodex Reversed Phase and GPC columns is known for strong chemical resistance, pH stability, and reliable performance characteristics, especially for polymer-based materials.

Interested? Contact us for current availability at info@shodex.de and mention the specific serial number when placing your order.

RP

Product code	Product name	Serial number	Qty	Functional group	Gel material	Plates per column	Particle size	Pore size	Size (ID x length)	Shipping solvent	USP	Info	Regular price	Sales price
F7620009	Asahipak ODP-50 2D	M1570039	1	C18	Polyvinyl alcohol	≥ 5,000	5 µm	250 Å	2.0 x 150 mm	H2O/ACN = 35/65	L67	Carbon load 17%	1.030,00 EUR	340,00 EUR
F6713001	Asahipak ODP-50G 2A	G34T0073	1	C18	Polyvinyl alcohol	(guard)	5 µm	-	2.0 x 10 mm	H2O/ACN = 35/65	L67	Carbon load 17%	542,00 EUR	176,00 EUR
F6713001	Asahipak ODP-50G 2A	G34T0074	1	C18	Polyvinyl alcohol	(guard)	5 µm	-	2.0 x 10 mm	H2O/ACN = 35/65	L67	Carbon load 17%	542,00 EUR	176,00 EUR
F6650040	Silica C18M 4D	K1242014	1	C18	Silica	≥ 10,000	5 µm	100 Å	4.6 x 150 mm	H2O/CH3OH = 30/70	L1	Carbon load 16%	670,00 EUR	268,00 EUR
F6650040	Silica C18M 4D	K1242015	1	C18	Silica	≥ 10,000	5 µm	100 Å	4.6 x 150 mm	H2O/CH3OH = 30/70	L1	Carbon load 16%	670,00 EUR	268,00 EUR
F6650040	Silica C18M 4D	K1242016	1	C18	Silica	≥ 10,000	5 µm	100 Å	4.6 x 150 mm	H2O/CH3OH = 30/70	L1	Carbon load 16%	670,00 EUR	268,00 EUR
F6650040	Silica C18M 4D	K1242017	1	C18	Silica	≥ 10,000	5 µm	100 Å	4.6 x 150 mm	H2O/CH3OH = 30/70	L1	Carbon load 16%	670,00 EUR	268,00 EUR

GPC

Product code	Product name	Serial number	Qty	Gel material	Plates per column	Particle size	Pore size	Size (ID x length)	Shipping solvent	USP	MW range	Exclusion limit	Regular price	Sales price
F6025050	GPC HK-405	G67J3023	1	Styrene divinylbenzene copolymer	≥ 7,000	3 µm	5,000 Å	4.6 x 150 mm	THF	L21	10,000 - 2,500,000	4,000,000	2.434,00 EUR	886,00 EUR
F6025050	GPC HK-405	G67J3024	1	Styrene divinylbenzene copolymer	≥ 7,000	3 µm	5,000 Å	4.6 x 150 mm	THF	L21	10,000 - 2,500,000	4,000,000	2.434,00 EUR	886,00 EUR
F6028114	GPC KF-402.5HQ	E2110036	1	Styrene divinylbenzene copolymer	≥ 25,000	3 µm	300 Å	4.6 x 250 mm	THF	L21	300 - 10,000	20,000	1.904,00 EUR	691,00 EUR
F6021043	GPC LF-404	E2580002	1	Styrene divinylbenzene copolymer	≥ 14,000	6 µm	3,000 Å	4.0 x 250 mm	THF	L21	300 - 2,000,000	2,000,000	2.005,00 EUR	722,00 EUR