

Shodex Amino HILIC NH2P-50 Series

The **EXACT** replacement for your Grace Prevail Carbohydrate ES Column



Alltech product number	Alltech product name	Shodex equivalent product number	Shodex equivalent product name	Functional group	Particle size	Pore Size	Stationary phase	I.D.	Length
35102	Prevail ES 4.6 x150mm	F7630002	Asahipak NH2P-50 4D	NH2	5µm	100 A	Polyvinylalcohol	4.6 mm	150 mm
35101	Prevail ES 4.6 x250mm	F7630001	Asahipak NH2P-50 4E	NH2	5µm	100 A	Polyvinylalcohol	4.6 mm	250 mm
35116	Prevail ES 2.0 x150mm	F7630006	Asahipak NH2P-50 2D	NH2	5µm	100 A	Polyvinylalcohol	2.0 mm	150 mm
-	Prevail ES 2.0 x250mm	F7630010	Asahipak NH2P-40 2E	NH2	4µm	100 A	Polyvinylalcohol	2.0 mm	250 mm

[Please also check out other Shodex Carbohydrate columns HILICpak VG-50, VC-50, VN-50 and VT-50](#)



www.shodex.de

Showa Denko Europe GmbH, Konrad-Zuse-Platz 3, 81829 Munich, Germany,
www.shodex.de; info@shodex.de;
Tel. +49 (0) 89 93996234