


Shodex™



HPLC Columns

Shodex Protein Kits

- Suitable for orthogonal analysis of proteins with three different separation techniques, among others SEC, RPC and IEC (strong cation exchanger)
- Kits include three columns: one column from PROTEIN KW-800 series, RSpak RP18-415 and IEC SP-825
- Three different Shodex Protein Kits are available:


Kit name	Small Protein Kit (up to 150 kDa)	Medium Protein Kit (up to 1000 kDa)	Large Protein Kit (up to 4000 kDa)
Article number	3240001	3240002	3240003
Columns	1. Protein KW-802.5 2. RSpak RP18-415 3. IEC SP-825	1. Protein KW-803 2. RSpak RP18-415 3. IEC SP-825	1. Protein KW-804 2. RSpak RP18-415 3. IEC SP-825

SHOWA
DENKO
EUROPE

Showa Denko Europe GmbH
Shodex Business
Konrad-Zuse-Platz 3
81829 Munich, Germany

+49 (0) 89 93 99 62 37
info@shodex.de
www.shodex.de
Europe, Middle East,
Africa, Russia

Shodex™



HPLC Columns

Shodex MALS Antibody Kit

- Suitable for high resolution analysis of antibodies (e.g. Rs (IgG Monomer/Dimer) = 2.2) and proteins
- Kit includes two low bleeding SEC columns: PROTEIN LW-803 and PROTEIN LW-403 4D
- SEC columns with target molecular weight range between 10 and 700 kDa (Exclusion Limit: 1000 kDa)
- PROTEIN LW-403 4D: sub-2 µm particles for fast analysis


Kit name	MALS Antibody Kit
Article number	3240004
Columns	1. Protein LW-803 2. Protein LW-403 4D

SHOWA
DENKO
EUROPE

Showa Denko Europe GmbH
Shodex Business
Konrad-Zuse-Platz 3
81829 Munich, Germany

+49 (0) 89 93 99 62 37
info@shodex.de
www.shodex.de
Europe, Middle East,
Africa, Russia

Shodex™



HPLC Columns

Shodex Method Development Kits

- Suitable for the screening and development of new methods for unknown samples with different polarities
- Kit includes three columns (4.6 x 250 mm) which cover following separation techniques: RPC and HILIC. It is also available in smaller column size, 2.0 x 150 mm, for LC/MS applications:

Kit name	Method Development Kit	Method Development Kit LC/MS
Article number	3240005	3240006
Columns	1. Asahipak ODP-50 4E 2. ODP2 HP-4E 3. Asahipak NH2P-50 4E	1. Asahipak ODP-50 2D 2. ODP2 HP-2D 3. Asahipak NH2P-50 2D

SHOWA
DENKO
EUROPE

Showa Denko Europe GmbH
Shodex Business
Konrad-Zuse-Platz 3
81829 Munich, Germany

+49 (0) 89 93 99 62 37
info@shodex.de
www.shodex.de
Europe, Middle East,
Africa, Russia